

Transparent addition-curing Silicone

# Fegura® Sil



Transparent silicone moulds made with **Fegura® Sil GLASS** are suitable for processing light-curing dental materials (e.g. composite resins). Light-curing through the silicone mould is advantageous.

The use of **Fegura® Sil GLASS** silicone moulds highly reduces the oxygen inhibit layer of composites. Thus the surface of the composite is denser and can be polished easily.

**Fegura® Sil GLASS** is available in two versions with different hardness:

- ✓ Shore A 55 **regular**
- ✓ Shore A 75 **super hard**

REF	PRODUCT	QUANTITY
FSG0120	Fegura® Sil Glass regular Shore A 55	2 x 50 ml cartridges 12 mixer yellow
FSG0100	Fegura® Sil Glass super hard Shore A 75	2 x 50 ml cartridges 12 mixer yellow



**Feguramed GmbH** • Jahnstraße 2 • D-74722 Buchen-Hettingen • GERMANY  
Tel. +49 (0)6281 / 5227-0 • Fax +49 (0)6281 / 5227-15  
E-Mail: [info@feguramed.com](mailto:info@feguramed.com) • Internet: [www.feguramed.com](http://www.feguramed.com)