

HIGH-END POLISHING PASTE

# Sunny Boy MO



PLAYING THE BEST - WITH DIAMOND

**High-End Polishing Paste for Non-Precious Alloys**

## Sunny Boy MO

An exceptionally dense, high-lustrous surface of dental alloys is of utmost importance. It is prerequisite for the avoidance of plaque deposits and thus for a longer lifetime of dentures. Feguramed's Sunny Boy MO is an abrasive polishing paste, particularly suitable for CoCr- and NiCr alloys. Perfect high lustre from pre-polishing to mirror finishing in a single working step is possible, because of the optimized composition of abrasives. Sunny Boy MO is wax-bound and ensures clean and almost dust-free processing. Thus the paste adheres extremely good to the polishing brush and offers economical consumption.

### SALES UNIT:

Bar: 550 g  
Color: green  
Ref.: 8073



**Feguramed GmbH** • Jahnstraße 2 • D-74722 Buchen-Hettingen • GERMANY  
Tel. +49 (0)6281 / 5227-0 • Fax +49 (0)6281 / 5227-15  
E-Mail: [info@feguramed.com](mailto:info@feguramed.com) • Internet: [www.feguramed.com](http://www.feguramed.com)