


The new duplicating silicone

Fegura® Sil



IF FAST DUPLICATING IS IN DEMAND...

... TIME IS MONEY!

The new duplicating silicone

Fegura® Sil quick

Offers a very short setting time and is twice as fast as conventional duplicating silicones! As is the case for all Fegura® Sil duplicating silicones we guarantee outstanding quality criteria concerning precision, flowability, dimensional stability and tear resistance.



Color: A = white
B = orange

Shore A 18-20°

Working time: approx. 3 min (22° C)

Demoulding: after approx. 12-15 min (22° C)

SALES UNITS	REF.
1 kg A + 1 kg B	1080
6 kg A + 6 kg B	1081
25 kg A + 25 kg B	1082



Feguramed GmbH • Jahnstraße 2 • D-74722 Buchen-Hettingen • GERMANY
Tel. +49 (0)6281 / 5227-0 • Fax +49 (0)6281 / 5227-15
E-Mail: info@feguramed.com • Internet: www.feguramed.com